


WATER FLOSS PRO

CORDLESS DENTAL WATER JET

IMPORTANT:

READ ALL INSTRUCTIONS BEFORE USE

Before use, you must prepare the Water Floss Pro. Remove the reservoir and place the Pump Housing, Tube and Nozzle completely under water. Pump the button for 60 seconds. This will allow the unit to function properly.



To Operate:

- 1. Fill the reservoir with water. Fill to the "Max line" on the reservoir. Do not over fill.**
- 2. Press Button to water floss.**