


2.



## Circulation Leg Exerciser Instructions

1. Connect power adapter to the unit and power jack to a power outlet.
2. Slide the power switch (B) to 1 for slow, 2 for normal speed, or to OFF.
3. Sit on a chair and place your feet on the pedals (A).
4. When finish, turn unit off.

### Optional

To angle the unit upward, pull the legs (C) from the underside.

### Caution

- Do not stand or sit on the unit.
- Keep out of reach of children. This is not a toy.
- Do not use on slippery or uneven surfaces.
- Do not wet the unit.
- Do not place your fingers or toes between the pedals.

1.

